

**Product Name:** MAGEA1**Product:** Purified Recombinant Protein**Lot Number:****Production Date****Catalog No.:** DP0061-P1**Concentration:****Expiring Date:****Protein Purification****Starting Material:** Cell Lysate**Buffer/Components:** 25 mM Tris, 500 mM NaCl, 0.4 mg/ml Triple Flag, pH 7.3, 10% glycerol**Storage Conditions:** Store at -80°C. Avoid repeated freeze-thaw cycles.

<b>Test:</b>	<b>Method:</b>	<b>Document:</b>	<b>Results:</b>
Protein Concentration	Bradford protein assay method	G3476	
Purity	SDS-Page (reduced)	G3302	

**Prepared By:** \_\_\_\_\_

Quality Control-Date

**Reviewed By:** \_\_\_\_\_

Quality Assurance-Date