

Product Name: SOX2
Product: Purified Recombinant Protein
Lot Number:
Production

Catalog No.: DP0062-P1
Concentration:
Expiring Date:

Protein Purification

Starting Material: Cell Lysate
Buffer/Components: 25 mM Tris-HCl, 100 mM glycine, 10% glycerol, pH 7.3
Storage Conditions: Store at -80°C

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method		
Purity	SDS-Page (reduced)	G3302	

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date