

**Product Name:** IL-6  
**Product:** Purified Recombinant Protein  
**Lot Number:**  
**Production Date:**

**Catalog No.:** DP0063-P1  
**Concentration:**  
**Expiring Date:**

**Protein Purification**

**Starting Material:** cell culture supernatant  
**Buffer/Components:** PBS, pH 7.4  
**Storage Conditions:** Store at -80°C. Avoid repeated freeze-thaw cycles.

<b>Test:</b>	<b>Method:</b>	<b>Document:</b>	<b>Results:</b>
Protein Concentration	Microplate BCA method	G3476	
Purity	SDS-Page (reduced)	G3302	

**Prepared By:** \_\_\_\_\_  
Quality Control-Date

**Reviewed By:** \_\_\_\_\_  
Quality Assurance-Date