

Product Name: CENP-B**Product:** Purified Recombinant Protein**Lot Number:****Production Date:****Catalog No.:** DP0066-P1**Concentration:****Expiring Date:****Protein Purification****Starting Material:** Cell Lysate**Buffer/Components:** 25mM Tris-HCl, 0.4mg/ml DDK, 500mM NaCl, 10% glycerol, pH8.5**Storage Conditions:** Store at -80°C. Avoid repeated freeze-thaw cycles.

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method		
Purity	SDS-Page (reduced)	G3302	

Prepared By: _____

Quality Control-Date

Reviewed By: _____

Quality Assurance-Date