

**Product Name:** GBU4-5**Product:** Purified Recombinant Protein**Lot Number:****Production Date:****Catalog No.:** DP0068-P1**Concentration:****Expiring Date:****Protein Purification****Starting Material:** Cell Lysate**Buffer/Components:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol**Storage Conditions:** Store at -80°C. Avoid repeated freeze-thaw cycles.

<b>Test:</b>	<b>Method:</b>	<b>Document:</b>	<b>Results:</b>
Protein Concentration	Bradford protein assay method	G3476	
Purity	SDS-Page (reduced)	G3302	> 80%

**Prepared By:** \_\_\_\_\_

Quality Control-Date

**Reviewed By:** \_\_\_\_\_

Quality Assurance-Date