

Product Name: TRIM33**Product:** Purified Recombinant Protein**Lot Number:****Production Date:****Catalog No.:** DP0069-P1**Concentration:****Expiring Date:****Protein Purification****Starting Material:** Cell Lysate**Buffer/Components:** 25 mM Tris-HCl, 100 mM glycine, 250 mM NaCl, 10% glycerol, pH 7.3**Storage Conditions:** Store at -80°C. Avoid repeated freeze-thaw cycles.

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method	G3476	
Purity	SDS-Page (reduced)	G3302	

Prepared By: _____

Quality Control-Date

Reviewed By: _____

Quality Assurance-Date