

**Product Name:** PCT  
**Product:** Purified Recombinant Protein      **Catalog No.:** TP750123  
**Lot Number:**      **Concentration :**  
**Manufactured:**      **EXP date:**

### Protein Purification

**Starting Material:** Lysate  
**Buffer/Components:** PBS, pH7.4, 10% glycerol.  
**Storage Conditions:** -80°C

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method		
Purity	SDS-Page (reduced)	G3302	

**Certified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_