


Lenti RapidTiter Sticks Reference Card



Cat #: EB200001/EB200002/EB200003

Lot #: AR239801

Optimal Detection Range



$> 2.5 \times 10^7$
TU/mL

2.5×10^7
TU/mL

1.0×10^7
TU/mL


2.5×10^6
TU/mL

1.0×10^6
TU/mL


2.5×10^5
TU/mL

$< 2.5 \times 10^5$
TU/mL


Estimate Your Titer in 4 Easy Steps




**Add 20ul of Your
Sample to
Sample Well**



**Add 80ul of
Running Buffer
to Sample Well**



**Wait for 10
mins**



**Read
Results**

Scan here for the
full product manual

