

Antibody Diluent Solution

Storage: 2-8°C

Catalog No: E09-500 500 ml
E09-300 60mlx5
E09-Sample 50ml

Intended Use:

Antibody Diluent Solution is used to dilute primary and secondary antibodies. The solution is not suitable for those antibodies conjugated to peroxidase enzyme because NaN_3 in the diluent will inhibit the activity of peroxidase.

Kit Components:

Catalog No.	Size
E09-500	500 ml
E09-300	60mlx5

Antibody Diluent Solution contains 0.05% NaN_3 .

Recommended Protocol:

Refer the dilution factors on manufacture's specification sheet or insert. For example if dilution is 50x, add 1ml of the concentrated antibody into 50ml of Antibody Diluent to obtain working dilution. Golden Bridge International's Antibody Diluent Solution may be used for antibodies from other source.

Precaution:

Wear gloves and take other necessary laboratory safety procedure.

Storage:

Store at 2-8°C.

Remarks:

For research use only.