

Product datasheet for 20061

Tacr3 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Reactivity:	Canine, Gerbil, Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Rattus norvegicus tachykinin receptor 3 (Tacr3)
Database Link:	Entrez Gene 6870 Human Entrez Gene 21338 Mouse Entrez Gene 24808 Rat
Synonyms:	TAC3R, Neuromedin-K receptor, NKR, Neurokinin B receptor, NK-3 receptor, NK-3R, NK3R

Product images:




Staining with Neurokinin 3 Receptor Rabbit Antibody (20061)




Staining with Neurokinin 3 Receptor Rabbit Antibody (20061)



[View online »](#)



Staining with Neurokinin 3 Receptor Rabbit Antibody (20061)



Western blot with Neurokinin 3 Receptor Rabbit Antibody (20061)