


## Product datasheet for 20064

### TAC1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens tachykinin precursor 1 (TAC1), transcript variant beta
Database Link:	<a href="#">Entrez Gene 6863 Human</a>
Synonyms:	Substance-P
Protein Families:	Druggable Genome, Secreted Protein

#### Product images:



Staining with Substance P antibody



[View online »](#)



Staining with Substance P antibody