


Product datasheet for 20068

Oxytocin neurophysin 1 (OXT) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic oxytocin coupled to bovine thyroglobulin (BTg) with carbodiimide (CDI) linker.
Specificity:	Reacts With: Bat, Bovine, Chicken, Dog, Frog (<i>Xenopus Laevis</i>), Guinea Pig, Hamster, Human, <i>Manduca Sexta</i> (Moth), <i>Melopsittacus Undulatus</i> (Bird), Mole Rat, Monkey, Mouse, Rat, Sheep, Vole
Formulation:	State: Serum State: Lyophilized Whole Serum Preservative: ≤ 0.09% sodium azide
Gene Name:	Homo sapiens oxytocin/neurophysin I prepropeptide (OXT)
Database Link:	Entrez Gene 5020 Human
Synonyms:	OXT, Ocytocin
Protein Families:	Secreted Protein

Product images:



Staining with Oxytocin Antibody 20068



[View online »](#)