


Product datasheet for 20069

Vasopressin (AVP) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Arginine Vasopressin coupled to bovine thyroglobulin with carbodiimide.
Specificity:	Chicken, Hamster, Mouse, Rat.
Formulation:	State: Serum State: Lyophilized Whole Serum Preservative: $\leq 0.09\%$ sodium azide
Gene Name:	arginine vasopressin
Synonyms:	Arginine vasopressin, Antidiuretic Hormone, AVP

Product images:



20069 Vasopressin antibody staining of Rat Brain.



[View online »](#)