

Product datasheet for 20070

ACTH (POMC) (1-39) Rabbit Polyclonal Antibody


Product data:

| | |
|-------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | IF, IHC |
| Reactivity: | Rat |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | Adrenocorticotrophic hormone (ACTH) corresponding to amino acids (1-39) generated against porcine pituitary |
| Specificity: | Reacts With: Dove, Eel, Goldfish, Hamster, Human, Monkey, Mouse, Ram (Sheep), Rat, Trout. |
| Formulation: | State: Serum State: Lyophilized whole serum Preservative: ≤ 0.09% sodium azide |
| Gene Name: | Homo sapiens proopiomelanocortin (POMC), transcript variant 2 |
| Database Link: | Entrez Gene 24664 Rat |
| Synonyms: | Corticotropin, Adrenocorticotrophic hormone, POMC |
| Protein Families: | Druggable Genome |
| Protein Pathways: | Adipocytokine signaling pathway, Melanogenesis |



[View online »](#)

Product images:



Staining with ACTH (Adrenocorticotrophic Hormone) Antibody (20070)