


## Product datasheet for 20077

### VIP Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Porcine VIP coupled to bovine thyroglobulin (BTg) with carbodiimide (CDI) linker.
Formulation:	State: Serum State: Lyophilized whole serum (100 $\mu$ L) Preservative: $\leq$ 0.09% sodium azide
Gene Name:	vasoactive intestinal peptide
Synonyms:	VIP, Vasoactive intestinal peptide

#### Product images:



Staining of rat cortex with VIP (Vasoactive Intestinal Peptide) Antibody 20077



[View online »](#)