


Product datasheet for 20085

Chromogranin A (CHGA) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens chromogranin A (CHGA), transcript variant 1
Database Link:	Entrez Gene 1113 Human
Synonyms:	Chromogranin-A, CgA, CHGA, SP-1
Protein Families:	Druggable Genome, Secreted Protein

Product images:



Staining with SP-1 Chromogranin A (Bovine)
Antibody 20085



[View online »](#)