


## Product datasheet for 20086

### Chromogranin A (CHGA) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens chromogranin A (CHGA), transcript variant 1
Database Link:	<a href="#">Entrez Gene 1113 Human</a>
Synonyms:	Chromogranin-A, CgA, CHGA, SP-1
Protein Families:	Druggable Genome, Secreted Protein

#### Product images:



Staining with SP-1 Chromogranin A (Porcine)  
Antibody 20086



[View online »](#)