


Product datasheet for 22806

Dopamine beta Hydroxylase (DBH) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Bovine, Equine, Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	dopamine beta-hydroxylase
Database Link:	Entrez Gene 1621 Human
Synonyms:	DBH, Dopamine beta-monoxygenase, DBM, DOPO

Product images:



Staining with DBH (Dopamine-beta-Hydroxylase) Antibody 22806



[View online »](#)