


Product datasheet for 24427

Calbindin (CALB1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Bovine, Human, Monkey, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	calbindin 1
Database Link:	Entrez Gene 793 Human
Synonyms:	CALB1, CAB27, Calbindin D28, D-28K

Product images:



Calbindin D-28 Antibody (24427)



[View online »](#)