


Product datasheet for 24429

5HT5A receptor (HTR5A) (17-34) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens 5-hydroxytryptamine receptor 5A (HTR5A)
Database Link:	Entrez Gene 15563 Mouse Entrez Gene 25689 Rat
Synonyms:	5-hydroxytryptamine receptor 5A, 5-HT5A, 5-HT-5A, 5-HT-5
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Calcium signaling pathway, Neuroactive ligand-receptor interaction

Product images:



Staining with 5-HT (Serotonin) 5A Receptor Antibody (24429)



[View online »](#)