


Product datasheet for 24430

5HT7 Receptor (HTR7) (8-23) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens 5-hydroxytryptamine receptor 7 (HTR7), transcript variant a
Database Link:	Entrez Gene 15566 MouseEntrez Gene 65032 RatEntrez Gene 3363 Human
Synonyms:	5-hydroxytryptamine receptor 7, 5-HT-7, 5HT7
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Calcium signaling pathway, Neuroactive ligand-receptor interaction

Product images:



Staining with 5-HT (Serotonin) 7 Receptor
Antibody 24430



[View online »](#)