

## Product datasheet for 24505

### HTR2C Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens 5-hydroxytryptamine (serotonin) receptor 2C, G protein-coupled (HTR2C), transcript variant 2
Database Link:	<a href="#">Entrez Gene 15560 Mouse</a> <a href="#">Entrez Gene 25187 Rat</a> <a href="#">Entrez Gene 3358 Human</a>
Synonyms:	HTR1C, 5HT-1C, 5-hydroxytryptamine receptor 2C, 5-HTR2C, 5-HT-2C, 5-HT2C
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Calcium signaling pathway, Gap junction, Neuroactive ligand-receptor interaction

#### Product images:



Staining with 5-HT (Serotonin) 2C Receptor Antibody (24505)



[View online »](#)