

Product datasheet for 24506


NPY1R (356-382) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Frog, Guinea Pig, Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	neuropeptide Y receptor Y1
Database Link:	Entrez Gene 4886 Human
Synonyms:	NPYR, NPYY1, Neuropeptide Y receptor 1

Product images:


Slide of Rat Cortex at 400x



Staining with Neuropeptide Y Y1 Receptor Antibody (24506)



[View online »](#)



Staining with Neuropeptide Y Y1 Receptor Antibody (24506)