

## Product datasheet for 2AN6-Adn27

### Adiponectin (ADIPOQ) Mouse Monoclonal Antibody [Clone ID: Adn27]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	Adn27
Applications:	ELISA, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	adiponectin, C1Q and collagen domain containing
Database Link:	<a href="#">Entrez Gene 9370 Human</a>
Synonyms:	ADIPOQ, ACDC, ACRP30, APM1, GBP28



[View online »](#)