

Product datasheet for 2AN6-Adn36

Adiponectin (ADIPOQ) Mouse Monoclonal Antibody [Clone ID: Adn36]

Product data:

Product Type:	Primary Antibodies
Clone Name:	Adn36
Applications:	ELISA, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	adiponectin, C1Q and collagen domain containing
Database Link:	Entrez Gene 9370 Human
Synonyms:	ADIPOQ, ACDC, ACRP30, APM1, GBP28



[View online »](#)