

Product datasheet for **2AN6-Adn63**

Adiponectin (ADIPOQ) Mouse Monoclonal Antibody [Clone ID: Adn63]

Product data:

Product Type:	Primary Antibodies
Clone Name:	Adn63
Applications:	ELISA, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens adiponectin, C1Q and collagen domain containing (ADIPOQ), transcript variant 2
Database Link:	Entrez Gene 9370 Human
Synonyms:	ADIPOQ, ACDC, ACRP30, APM1, GBP28
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Adipocytokine signaling pathway, PPAR signaling pathway, Type II diabetes mellitus



[View online »](#)