

Product datasheet for 2G2-GhB9

Growth Hormone (GH1) Mouse Monoclonal Antibody [Clone ID: GhB9]

Product data:

Product Type:	Primary Antibodies
Clone Name:	GhB9
Applications:	ELISA
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens growth hormone 1 (GH1), transcript variant 1
Database Link:	Entrez Gene 2688 Human
Synonyms:	Pituitary growth hormone, Growth hormone 1, HGH
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway, Neuroactive ligand-receptor interaction



[View online »](#)