

Product datasheet for 2H8-F1

hCG beta (CGB3) (alpha chain) Mouse Monoclonal Antibody [Clone ID: F1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	F1
Reactivity:	Human
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens chorionic gonadotropin beta subunit 3 (CGB3)
Database Link:	Entrez Gene 1082 Human
Synonyms:	Choriogonadotropin subunit beta, CGB, CGB3, beta hCG, hCG-beta
Protein Families:	Druggable Genome



[View online »](#)