

Product datasheet for 2I1-3A6

Insulin (INS) (+Proinsulin) Mouse Monoclonal Antibody [Clone ID: 3A6]

Product data:

Product Type:	Primary Antibodies
Clone Name:	3A6
Applications:	ELISA
Reactivity:	Bovine, Human, Porcine
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens insulin (INS), transcript variant 1
Database Link:	Entrez Gene 3630 Human
Synonyms:	INS
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein
Protein Pathways:	Insulin signaling pathway, Maturity onset diabetes of the young, mTOR signaling pathway, Oocyte meiosis, Progesterone-mediated oocyte maturation, Prostate cancer, Regulation of actin cytoskeleton, Regulation of autophagy, Type I diabetes mellitus, Type II diabetes mellitus



[View online »](#)