

Product datasheet for 2LH2-L1

Luteinizing Hormone beta (LHB) Mouse Monoclonal Antibody [Clone ID: L1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	L1
Applications:	ELISA, WB
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens luteinizing hormone beta polypeptide (LHB)
Database Link:	Entrez Gene 3972 Human
Synonyms:	LHB, LH-B, LSH-beta, Luteinizing hormone subunit beta, Lutropin beta chain, Lutropin subunit beta
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	GnRH signaling pathway, Neuroactive ligand-receptor interaction



[View online »](#)