

Product datasheet for **4C28-CRP169**

C Reactive Protein (CRP) Mouse Monoclonal Antibody [Clone ID: CRP169]

Product data:

| | |
|----------------|---|
| Product Type: | Primary Antibodies |
| Clone Name: | CRP169 |
| Applications: | ELISA, IF, WB |
| Reactivity: | Human |
| Host: | Mouse |
| Isotype: | IgG2a |
| Clonality: | Monoclonal |
| Formulation: | State: Purified |
| Gene Name: | C-reactive protein, pentraxin-related |
| Database Link: | Entrez Gene 1401 Human |
| Synonyms: | PTX1, C Reactive Protein, Pentraxin-related |



[View online »](#)