

## Product datasheet for 4ER1-Epo2

### EPO Mouse Monoclonal Antibody [Clone ID: Epo2]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	Epo2
Applications:	ELISA
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens erythropoietin (EPO)
Database Link:	<a href="#">Entrez Gene 2056 Human</a>
Synonyms:	Epoetin
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway



[View online »](#)