

Product datasheet for 4F29-25

FABP3 Mouse Monoclonal Antibody [Clone ID: 25]

Product data:

Product Type:	Primary Antibodies
Clone Name:	25
Applications:	ELISA
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens fatty acid binding protein 3, muscle and heart (FABP3)
Database Link:	Entrez Gene 2170 Human
Synonyms:	FABP11, MDGI, H-FABP, Heart-type fatty acid-binding protein
Protein Pathways:	PPAR signaling pathway



[View online »](#)