

## Product datasheet for 4F29-30

### FABP3 Mouse Monoclonal Antibody [Clone ID: 30]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	30
Applications:	ELISA, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens fatty acid binding protein 3, muscle and heart (FABP3)
Database Link:	<a href="#">Entrez Gene 2170 Human</a>
Synonyms:	FABP11, MDGI, H-FABP, Heart-type fatty acid-binding protein
Protein Pathways:	PPAR signaling pathway



[View online »](#)