

Product datasheet for **4FBN3-FND5**

Fibronectin (FN1) Mouse Monoclonal Antibody [Clone ID: FND5]

Product data:

Product Type:	Primary Antibodies
Clone Name:	FND5
Applications:	ELISA, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens fibronectin 1 (FN1), transcript variant 3
Database Link:	Entrez Gene 2335 Human
Synonyms:	FN1, Cold-insoluble globulin, CIG
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein
Protein Pathways:	ECM-receptor interaction, Focal adhesion, Pathways in cancer, Regulation of actin cytoskeleton, Small cell lung cancer



[View online »](#)