

Product datasheet for 4GP31-1G6

GPBB (PYGB) Mouse Monoclonal Antibody [Clone ID: 1G6]

Product data:

Product Type:	Primary Antibodies
Clone Name:	1G6
Applications:	ELISA, WB
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens phosphorylase, glycogen; brain (PYGB)
Database Link:	Entrez Gene 5834 Human
Synonyms:	Glycogen phosphorylase brain form
Protein Families:	Druggable Genome
Protein Pathways:	Insulin signaling pathway, Starch and sucrose metabolism



[View online »](#)