

## Product datasheet for 4I22-GF1

### Interferon gamma (IFNG) Mouse Monoclonal Antibody [Clone ID: GF1]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	GF1
Applications:	ELISA, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens interferon, gamma (IFNG)
Database Link:	<a href="#">Entrez Gene 3458 Human</a>
Synonyms:	IFN-gamma, gamma IFN
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Allograft rejection, Cytokine-cytokine receptor interaction, Graft-versus-host disease, Jak-STAT signaling pathway, Natural killer cell mediated cytotoxicity, Proteasome, Regulation of autophagy, Systemic lupus erythematosus, T cell receptor signaling pathway, TGF-beta signaling pathway, Type I diabetes mellitus



[View online »](#)