

Product datasheet for 4IG8-G5F8

IGFBP1 Mouse Monoclonal Antibody [Clone ID: G5F8]

Product data:

Product Type:	Primary Antibodies
Clone Name:	G5F8
Applications:	ELISA, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens insulin like growth factor binding protein 1 (IGFBP1)
Database Link:	Entrez Gene 3484 Human
Synonyms:	IGF-binding protein 1, IGFBP-1, IBP1, IBP-1, Placental protein 12, PP12
Protein Families:	Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Secreted Protein



[View online »](#)