

Product datasheet for **4IGF4-IBP190**

IGFBP4 Mouse Monoclonal Antibody [Clone ID: IBP190]

Product data:

Product Type:	Primary Antibodies
Clone Name:	IBP190
Applications:	ELISA, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens insulin like growth factor binding protein 4 (IGFBP4)
Database Link:	Entrez Gene 3487 Human
Synonyms:	IGF-binding protein 4, IGFBP-4, IBP4, IBP-4
Protein Families:	Secreted Protein



[View online »](#)