

Product datasheet for **4LOX1-LOX19-22**

LOX 1 (OLR1) (sLOX-1) Mouse Monoclonal Antibody [Clone ID: LOX19-22]

Product data:

Product Type:	Primary Antibodies
Clone Name:	LOX19-22
Applications:	ELISA, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens oxidized low density lipoprotein (lectin-like) receptor 1 (OLR1), transcript variant 1
Database Link:	Entrez Gene 4973 Human
Synonyms:	LOX-1, Ox-LDL receptor 1, Ox-LDL receptor 1 soluble form, sLOX1, sLOX-1, CLEC8A
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
Protein Pathways:	PPAR signaling pathway



[View online »](#)