

Product datasheet for **4M40-M2F1**

EMA (MUC1) Mouse Monoclonal Antibody [Clone ID: M2F1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	M2F1
Applications:	ELISA, IF, IHC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	mucin 1, cell surface associated
Database Link:	Entrez Gene 4582 Human
Synonyms:	MUC-1, PEMT, Episialin, EMA, H23AG, PUM, DF3, CA 15-3



[View online »](#)