

Product datasheet for 4N6-5G10

NSE (ENO2) Mouse Monoclonal Antibody [Clone ID: 5G10]

Product data:

Product Type:	Primary Antibodies
Clone Name:	5G10
Applications:	ELISA, IHC, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	enolase 2
Database Link:	Entrez Gene 2026 Human
Synonyms:	NSE, ENO2, Enolase 2, Neural enolase, Gamma-enolase



[View online »](#)