

## Product datasheet for **4P41-10H9**

### **PAPP A (PAPPA) Mouse Monoclonal Antibody [Clone ID: 10H9]**

#### **Product data:**

|                |  |
|----------------|--|
| Product Type:  | Primary Antibodies                     |
| Clone Name:    | 10H9                                   |
| Applications:  | ELISA, IP, WB                          |
| Reactivity:    | Human                                  |
| Host:          | Mouse                                  |
| Isotype:       | IgG2a                                  |
| Clonality:     | Monoclonal                             |
| Formulation:   | State: Purified                        |
| Gene Name:     | pappalysin 1                           |
| Database Link: | <a href="#">Entrez Gene 5069 Human</a> |
| Synonyms:      | PAPP-A, IGFBP-4ase                     |



[View online »](#)