

## Product datasheet for **4P41-3C8**

### **PAPP A (PAPPA) Mouse Monoclonal Antibody [Clone ID: 3C8]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	3C8
Applications:	ELISA, IP, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	pappalysin 1
Database Link:	<a href="#">Entrez Gene 5069 Human</a>
Synonyms:	PAPP-A, IGFBP-4ase



[View online »](#)