

Product datasheet for **4PC47-44D9**

Calcitonin (CALCA) Mouse Monoclonal Antibody [Clone ID: 44D9]

Product data:

Product Type:	Primary Antibodies
Clone Name:	44D9
Applications:	ELISA, IF, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	calcitonin related polypeptide alpha
Database Link:	Entrez Gene 796 Human
Synonyms:	CALCA, CALC1



[View online »](#)