

Product datasheet for 4PD4-PAPP30

PAPP A (PAPPA) Mouse Monoclonal Antibody [Clone ID: PAPP30]

Product data:

Product Type:	Primary Antibodies
Clone Name:	PAPP30
Applications:	ELISA
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens pregnancy-associated plasma protein A, pappalysin 1 (PAPPA)
Database Link:	Entrez Gene 5069 Human
Synonyms:	PAPP-A, IGFBP-4ase
Protein Families:	Druggable Genome, Protease, Secreted Protein



[View online »](#)