

Product datasheet for **4RB2-RB45**

RBP4 Mouse Monoclonal Antibody [Clone ID: RB45]

Product data:

| | |
|-------------------|-----------------------------------------------------------------------------|
| Product Type: | Primary Antibodies |
| Clone Name: | RB45 |
| Applications: | ELISA |
| Reactivity: | Human |
| Host: | Mouse |
| Isotype: | IgG1 |
| Clonality: | Monoclonal |
| Formulation: | State: Purified |
| Gene Name: | Homo sapiens retinol binding protein 4 (RBP4), transcript variant 1 |
| Database Link: | Entrez Gene 5950 Human |
| Synonyms: | RBP-4, Retinol-binding protein 4, Plasma retinol-binding protein, PRBP, RBP |
| Protein Families: | Druggable Genome, Secreted Protein |



[View online »](#)