

Product datasheet for 4S37-3B10

S100A1 (ab) Mouse Monoclonal Antibody [Clone ID: 3B10]

Product data:

Product Type:	Primary Antibodies
Clone Name:	3B10
Applications:	ELISA, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens S100 calcium binding protein A1 (S100A1)
Database Link:	Entrez Gene 6271 Human
Synonyms:	S100-A1, S100A, S100 alpha, S100 calcium-binding protein A1, S-100 protein subunit alpha, S-100 protein alpha chain



[View online »](#)