

Product datasheet for 4T10-2C8

TNF alpha (TNF) Mouse Monoclonal Antibody [Clone ID: 2C8]

Product data:

Product Type:	Primary Antibodies
Clone Name:	2C8
Applications:	ELISA, FN, IHC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	tumor necrosis factor
Database Link:	Entrez Gene 7124 Human
Synonyms:	TNF, TNF-a, TNFA, TNFSF2, Cachectin



[View online »](#)