

Product datasheet for 4T15-1C10

Transferrin (TF) Mouse Monoclonal Antibody [Clone ID: 1C10]

Product data:

Product Type:	Primary Antibodies
Clone Name:	1C10
Applications:	ELISA, FN, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	transferrin
Database Link:	Entrez Gene 7018 Human
Synonyms:	Serotransferrin, Siderophilin



[View online »](#)